

<i>Type, chosen from dropdown menu (usually "Sample")</i>	<i>Autosampler location</i>	<i>Description</i>	<i>Must match filename of method that is loaded. Include everything but ".m".</i>	<i>The data file name. Must be unique or you will overwrite data.</i>
Type	Vial	Sample	Method/keyword	Data file
Sample	2	a-Pinene	chem472-18_method	group1_pinene001
Sample	2	a-Pinene replicate	chem472-18_method	group1_pinene002
Sample	5	Thujone	chem472-18_method	group1_thujone001
Sample	6	Limonene	chem472-18_method	group1_limonene001